BioSonic US200

Ultrasonic Scaler

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- EN** Owner's Guide
- FR** Guide d'utilisation
- ES** Guía de uso

Table of Contents

EN	Owner's Guide.....	4
FR	Guide d'utilisation	19
ES	Guía de uso.....	34

BioSonic US200 Ultrasonic Scaler

CAUTION: Federal law (United States) restricts this device to sale by or on the order of a licensed healthcare practitioner. This product is for use by trained healthcare professionals only.

 WARNING	To reduce the risk of patient or personal injury or delay of procedure, read this manual in its entirety before installing and operating this device. Improper use may lead to patient or user injury.
---	---

The BioSonic US200 Ultrasonic Scaler is an enhanced third-generation unit with an intuitive design and easy-to-use interface which allows for precision, flexibility, and comfort.

Clinical Benefits of Scaling

Scaling is the process by which calculus is removed from the tooth surface. This can be accomplished using either hand instruments (such as hand scalers or curettes) or ultrasonic scalers. When using ultrasonic scalers, high frequency waves generated within the ultrasonic handpiece cause the tip of the insert to vibrate at 25,000 or 30,000 cycles per second. This action can remove calculus from tooth surfaces more effectively and more rapidly than hand instruments.

Intended Purpose

Indications for Use

Ultrasonic Applications for:

- All sub-gingival and supra-gingival scaling and routine prophylaxis procedures.
- Periodontal therapy

Contraindications

The BioSonic US200 Ultrasonic Scaler should not be used for/on:

- Restorative dental procedures involving the condensation of amalgam
- Surgical or endodontic procedures
- Patients highly susceptible to infections
- Patients with chronic respiratory problems
- Patients with dysphagia or swallowing difficulty with water flow
- Patients with cardiovascular disease with secondary pulmonary disease or chronic pulmonary disease due to a risk of aspiration of microorganisms (MOs) in biofilm

Intended Patient Population

The BioSonic US200 is intended for use on patients with primary to geriatric teeth requiring ultrasonic scaling.

I. Unpacking

A. Remove the unit and components from the shipping carton

The shipping carton should contain:

- [1] BioSonic US200 Unit
- [1] Handpiece Assembly
- [2] Spare Autoclavable Handpiece Sheaths
- [1] Footswitch Assembly
- [1] Water Line Assembly with Filter and Quick Disconnect Fitting
- [2] Spare Water Line Filters
- [2] Spare Handpiece O-Rings
- [1] Power Cord
- [1] Power Adapter
- [1] Instructions for Use Guide

- B. Inspect components for any damage which may have occurred during shipping.
- C. Check the serial number on the underside of the unit against the number on the shipping carton. If these numbers are not the same, contact the dealer.
- D. Complete the warranty within 10 days of receipt online at warranty.coltene.com.

II. Set Up

A. Placement

- Position the unit for easy access during patient treatment near a wall receptacle and cold-water line while ensuring there is enough space for the device to easily connect and disconnect from the outlet and water lines. Ensure the unit is placed on a solid flat surface, such as a countertop.

B. Water line connectors

- Incoming cold-water temperature should be between 5° C and 25° C (41° - 77° F). Supply line pressure to the unit must be between a minimum of 25 psi (172 kPa/1.724 bar) and a maximum of 50 psi (345 kPa/3.447 bar).
- The use of an in-line filter to be placed between the main shut-off valve and the unit is recommended. The main shut-off valve is a part of the plumbing in the operatory.
- Although the unit has a filter, it is not intended to replace a filter in the main line. The filter will capture smaller particles to further prevent handpiece clogging.
- The dental water line should be thoroughly flushed prior to connection to the unit.
- Only use waterline accessory provided with the unit or a replacement waterline provided by Dent4You AG.
- Always utilize the water inline filter provided with the unit.
- Dent4You AG recommends the use of water that meets EPA regulatory standards for drinking water (500 CFU/mL) for non-surgical routine dental procedures. Before use of the BioSonic US200 for the first time, an initial shock treatment should be conducted and followed by water testing to verify acceptable bacteria levels. These products should be used per specific manufacturer's instructions.

C. Water supply line connection

- Insert open end of the water line tube into the push-in water connector located on the rear of the unit. When pushing the water line into the push-in connector, there will be a bit of resistance. By pushing a little harder the line tubing will go past the internal O-ring making a complete seal. Pull back on the water tubing to ensure it is secure. (Fig. 1 & 2). If necessary, to disconnect the water line, hold the outer ring of the push-in water connector against the unit and pull the tube from the fitting.
- Check to ensure in-line filter is tightly attached to both sections of the tube.
- Connect end of tube with the quick-connect fitting to the office cold water supply line.
- Open the main shut-off valve.
- Inspect the unit and all connections to make certain there are no leaks, such as pooled water or dripping, prior to turning device on and beginning procedure.

D. Handpiece

- Push the connector on the handpiece assembly into the receptacle in the front of the unit. The connector will "click" into place when it is fully connected. The arrow on the connector indicates which side should be up when inserting the connector. (Fig. 3)
- Ensure the handpiece is securely connected.

E. Footswitch Assembly

- Push the footswitch connector into the receptacle in the back of the unit. The arrow on the connector indicates which side should be up when inserting the connector. (Fig. 4)

F. Power cord/power connection

- Plug the detachable power cord into the power adapter block. (Fig. 5)
- Plug the power adapter cord into the receptacle on the back of the unit. (Fig. 6)
- Insert the power plug of the detachable power cord into a GFCI outlet, preferably "Hospital Grade." To avoid risk of electrical shock, this equipment should only be connected to a supply mains with electrical earth.
- NOTE: If the power plug supplied with the unit does not fit the electric receptacle in the office, contact the dealer or an electrician.
- Check the LED below the power knob. If it is illuminated green, the unit is ON. If it is not illuminated, rotate the power knob clockwise to turn the unit on. The unit is shipped with the power knob in the OFF position.

III. Operating Instructions

A. Patient preparation and positioning

 WARNING	The scaling process carries the potential risk for aerosol transmission of infection. Be sure to follow appropriate infection control procedures as recommended by health authorities. Please refer to Section IV of this guide for Infection Control guidelines.
---	--

- Have the patient perform an antimicrobial pre-procedural mouth rinse to reduce bacterial contamination of the aerosol.
- Place the dental chair in a supine position for optimal access to both the maxillary and mandibular arches. This improves patient comfort and increases clinician visibility.
- Position the patient's head for optimum access to the quadrant and tooth surface being treated. Evacuate water and debris with a saliva ejector or high-volume evacuator.

B. Ultrasonic scaling procedures

NOTE: Refer to section IV (Maintenance and Infection Control) of these Instructions For Use for general procedures to follow daily and between patients.

General operating procedures:

1. Install the sterilized handpiece sheath by threading it onto the handpiece assembly until the tapered end is flush with the end of the handpiece. (Fig. 7) Scaler sheaths should be sterilized before first use on each patient per Section IV of this manual.
2. Be sure the main dental office water supply valve is open.
3. Turn the unit ON by rotating the power knob clockwise until a click is heard. Verify that the green LED below the power knob is lit. Note: Ensure that the footswitch is not depressed when

the unit is turned ON.

4. Turn the water control knob clockwise to the fully open position (3 full rotations).
5. Hold the handpiece (without an insert installed) over a sink or drain. Depress the footswitch and flush the water line for at least 2 minutes.
6. Select the sterilized scaler insert. Scaler inserts should be sterilized before first use on each patient as stated in Section IV. Inspect the insert for any wear or damage to the laminate stack or insert tip. Do not use a damaged insert.
7. Hold the handpiece upright over a sink or drain and press the footswitch until the handpiece is completely filled with water. This lubricates the insert O-ring and helps to ensure maximum cooling of the insert during operation.
8. Push an insert into the handpiece using a twisting motion to completely seat the insert into place. (Fig. 8)
9. Set the power control knob to the preferred operating positions. Power settings are numbered 1 through 10 – each increment represents a 10% change in overall power in "Normal" mode. NOTE: Use the lowest effective power setting possible when scaling to minimize patient discomfort and the risk of tooth surface defects.
10. Press the two-stage footswitch to the first position (for normal mode) or the second position (for turbo mode) and adjust the water control so that there is a fine mist of water spray.

NOTE: Prolonged use of the unit at high power or in turbo mode may require an increase in water flow to prevent overheating.

C. Turbo Mode

- The Turbo mode can be activated using the two-position footswitch in the second (fully depressed) position. (Fig. 9)
- When the Turbo mode is activated, the LED on the unit will change from green to orange.
- When the Turbo mode is activated, the output power will instantaneously increase to 150% of the current power setting or to the maximum power setting, whichever is lower. See the power dial for equivalent power setting when transitioning from normal to Turbo.
- To exit Turbo mode, release the footswitch to first position. The LED will return to green.
- The Turbo feature is helpful to remove stubborn calculus and is intended only to be used as a temporary mode of operation.
- It may be necessary to adjust the water control knob with the unit in Turbo mode to ensure adequate water delivery to the tip and the tooth. Prior to use on each patient, test the combination of intended power and water to ensure a proper water flow is achieved at the selected power setting.

D. Helpful hints

- Use the lowest effective power setting possible when scaling to minimize patient discomfort.
- If the handpiece becomes hot to the touch, there is either too little water, too much power or both. Adjust power and/or water settings accordingly.
- The rubber O-ring on the insert must be lubricated with water before placing it in the handpiece. Not lubricating the O-ring will make proper seating more difficult. Fully seat insert by pushing inward while twisting until seated. Do not apply

excessive force to insert. The grip on the insert should be flush against the handpiece.

- Check and verify that water is reaching the working end of the insert tip by holding the handpiece over a sink or drain and activating the unit. Adjust the water control knob to ensure adequate flow is achieved for the selected power setting. A fine mist of water should be coming from the scaler insert tip (Fig.10).
- If water fails to flow properly, check the supply line filter for clogging. Filters are disposable. Replace in case of discoloration, visible debris in filter, or reduced water flow. Do not attempt to unclog used filters.
- To prevent accidental contact with the lips, tongue, and cheek, the scaler insert should be activated just before contacting the tooth, but after the insert has been placed in the patient's mouth.
- It is best to hold the insert in a pen grasp, using a very light touch when scaling both supra and sub-gingivally. Consistent pressure should be maintained while scaling - too much pressure of the tip against the tooth surface will decrease the effectiveness of the insert, as tip vibration is obstructed. If more scaling power is required, the power level of the unit should be adjusted accordingly.
- Check the scaler insert and O-rings for wear prior to use on each patient.

- Develop an operating sequence in scaling to ensure thoroughness and to minimize adjustments in patient and/or operator position.

E. Patient discomfort can be caused by:

- Excessive pressure. Apply the side of the scaler tip gently to the tooth surface using minimal pressure.
- Incorrect tip positioning. Avoid contact of the scaler insert tip and the tooth. (Fig. 11) Direct the tip away from root surfaces and use side of scaler insert only. (Fig. 12)
- Not keeping the tip in motion on tooth. The tip should always be in motion during scaling. Use vertical, horizontal, or oblique overlapping strokes to increase scaling efficiency and decrease discomfort.
- Excessive power. Use of scaler in Turbo Mode for an extended period may lead to patient discomfort and tooth surface defects.
- If sensitivity persists, alternate treatment of sensitive areas and less sensitive areas and/or decrease power setting where possible.

F. Turn the unit OFF by rotating the power knob counterclockwise until a click is heard.

IV. Maintenance and Infection Control

Manufacturer: Dent4You AG **Method:** Low-level disinfection, steam sterilization **Symbol:** N/A

Devices: BioSonic US200 Base Unit, Ultrasonic Inserts, Handpiece Sheath.

These instructions only apply to BioSonic ultrasonic inserts manufactured by Dent4You AG

. For all other inserts, always refer to the manufacturer's instructions regarding maintenance and infection control procedures.

WARNINGS	Do not use cold sterilization, chemical disinfection, or other methods of sterilization, including chemical vapor and dry heat since these methods may severely damage the insert or sheath and will void the warranty. Do not use chemical disinfectants prior to sterilization or rapid deterioration of the materials may occur.
Limitations on processing	Inspect base unit and cords. If any damage is observed do not continue use, and replace items instead. Follow instructions under Section VII. Order Information Follow insert manufacturer's instructions regarding tip wear and discontinue further reprocessing and use of any inserts that do not meet the minimum requirements for tip wear. Inspect sheath for cracks or other damage. If any damage is found, discard and use a new sheath.

Instructions	
Initial treatment at the point of use	At the end of a procedure, flush the ultrasonic insert by leaving the insert in the handpiece and turning the water flow knob to produce a fine mist. Activate the water flow by depressing the footswitch, as done during normal operation, for a minimum of 1 minute. No additional treatment required at point of use.
Preparation before cleaning	Remove and rinse insert and sheath under running cold utility water to remove all visible debris.

Instructions		
Cleaning	Manual	Fully immerse the insert(s) and sheath(s) in a BioSonic General Purpose or Enzymatic Cleaner detergent solution prepared per the manufacturer's recommendation. Soak the components in the detergent solution for 5 minutes. Brush the insert(s) and sheath(s) thoroughly with a soft bristled brush. Pay close attention to all crevices of the device. Ensure all visible soil is removed. Rinse the insert(s) and sheath(s) under running warm utility water until all visual evidence of detergent is removed. Pay close attention to all cannulations and small crevices of the device when rinsing, and if necessary rinse these areas with a syringe or other small instrument to ensure all remaining solution is removed from the surface of the part.
Cleaning	Automated Ultrasonic Cleaner	Fully immerse the insert(s) and sheath(s) in a BioSonic General Purpose Cleaner or Enzymatic detergent solution prepared per the manufacturer's recommendation. Soak the components in the detergent solution for 5 minutes. Brush the insert(s) and sheath thoroughly with a soft bristled brush. Pay close attention to all crevices of the device. Ensure all visible soil is removed. Clean the insert(s) and sheath(s) in an ultrasonic cleaner (such as COLTENE's BioSonic Ultrasonic Cleaner) with BioSonic UC30 General Purpose Cleaner diluted in accordance with manufacturer's directions. Allow the devices to remain in the ultrasonic cleaner for a minimum of 5 minutes. Rinse the insert(s) and sheath(s) under running warm utility water until all visual evidence of detergent is removed. Pay close attention to all cannulations and small crevices of the device when rinsing and if necessary rinse these areas with a syringe or other small instrument to ensure all remaining solution is removed from the surface of the part.
Disinfection		Wipe base unit and handpiece cord with OPTIM 1 wipes. No additional disinfection of inserts or handpiece sheaths required prior to sterilization. DO NOT use chemical disinfectants prior to sterilization or rapid deterioration of the materials may occur.
Drying		Allow disinfected surfaces to air dry. Manually wipe dry inserts and sheaths with a clean, lint-free cloth.
Maintenance, Inspection and Testing		Periodically check and inspect the power cord, water lines, filter, footswitch, handpiece, and handpiece O-ring, and insert O-ring for damage, wear or clogging. Replace if necessary. Check connectors at both ends of the cord to make sure they are fully seated. Dent4You AG recommends that the waterline system be cleaned and disinfected by performing a shock treatment utilizing FDA cleared dental unit waterline cleaner and dispensing system at least weekly followed by a water test to verify acceptable bacteria levels. These products should be used per specific manufacturer's instructions. Inspect rubber O-ring on handpiece for wear and/or damage. When the O-ring is worn and/or damaged, leakage may occur. Inspect devices for visible soil. If visible soil remains, repeat the cleaning process. Inspect BioSonic insert O-ring for wear. Replace if necessary. Inspect sheath for cracks or other damage. If any damage is found, discard and use a new sheath.

Instructions	
Sterilization	Do not sterilize the base unit or accessories. Only the handpiece sheath and ultrasonic inserts are intended to be sterilized. Pre-Vacuum: Place cleaned and dried BioSonic Inserts and Sheaths in an FDA-cleared sterilizer pouch. Sterilize the pouched BioSonic inserts and sheaths via a pre-vac steam sterilizing cycle at 132°C for 4 minutes with a 30-minute drying time, or alternatively 134°C for 3 minutes with a 30-minute drying time. Do not use cold sterilization, chemical disinfection, or other methods of sterilization, including chemical vapor and dry heat, since these methods may severely damage the insert and will void the warranty.
Manufacturer Contact	service@coltene.com

Date issued: June 2024

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.


V. Transportation/Environmental/Storage Conditions

Environmental Conditions


• Intended for indoor use

• Maximum altitude: 2000m (6600 ft.)

• Temperature range:

• Relative humidity:


• Duty cycle: 4min. ON / 1.5min OFF


Storage Conditions

The unit must be stored and transported within the following conditions:


• Temperature range:

• Relative humidity:

• Atmospheric pressure:

VI. Information for Safety

 WARNING!	• Use of accessories provided by Dent4You AG are intended for use with this equipment. Any modifications may result in the BioSonic US200 Ultrasonic Scaler malfunctioning or cause a hazardous situation. However, ultrasonic inserts (both 25 kHz and 30 kHz) manufactured by Cavitron and Hu-Friedy may also be used. • Failure to inspect prior to use as instructed per this manual may lead to risk of patient or user injury, and may result in loss of functionality and delay of treatment. • Inspect unit for any evidence of water leaks, such as pooled water and dripping, prior to turning device on and beginning procedure. Failure to inspect prior to use may result in loss of functionality and delay of treatment. • Device is not intended for use in an oxygen-rich environment. Use in an oxygen-rich environment may lead to an increased risk of fire, causing burns.
--	---

 WARNING!	• All components not explicitly stated as being repairable and replaceable by the end-user require that the device be sent back to the manufacturer for repair or replacement. Any unauthorized modifications may result in the BioSonic US200 Ultrasonic Scaler malfunctioning or cause a hazardous situation. • The scaling process carries the potential risk for aerosol transmission of infection. Be sure to follow appropriate infection control procedures as recommended by health authorities. Please refer to Section IV of this guide for Infection Control guidelines. RISK OF LACERATION: • Avoid contact of the scaler insert tip and the tooth or soft tissue. Direct the tip away from root surfaces and use side of scaler insert only. Failure to properly place insert prior to activation may lead to an increased risk of laceration. • Periodically check and inspect the handpiece and base unit for damage or wear. If excessive wear is evident, or the insert has been bent, discard the insert immediately. It may break during use. • The scaling process carries the potential risk for aerosol transmission of infection. Be sure to follow appropriate infection control procedures as recommended by health authorities. Please refer to Section IV of this guide for Infection Control guidelines. • Use retraction to protect the tongue, cheek and lips when the scaler is in use to avoid contact between the two. RISK OF ELECTROMAGNETIC INTERFERENCE: • Radiation interferes with operation of implanted medical devices such as pacemakers, causing chest pain and shortness of breath. • Use of this device adjacent to or stacked with other equipment should be avoided because it can result in improper operation. • Use of accessories and cables other than those provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.
--	--

 WARNING!	To reduce the risk of minor burns or temperature-related discomfort to the patient or user: • Do not operate the unit without water or with too little water. • Use cold (5-25°C) (41 -77°F) water only. • Do not touch the metal working tip of an activated insert with bare hands. • Do not operate the device for more than 4 minutes on and less than 1.5 minutes off
 WARNING!	To reduce the risk of electric shock or loss of functionality: • Do not modify the power cord. • Do not immerse the unit in water. • Use of a Ground Fault Circuit Interrupt electrical outlet is recommended.
 WARNING!	To reduce the risk of water-borne illness: • Flush water supply to this device regularly. To reduce the risk of cross-contamination: Remove the sterilizable handpiece sheath and all ultrasonic inserts after use. Clean and sterilize between each patient, according to the instructions for use.
 CAUTION!	• Persons fitted with cardiac pacemakers or any other intra-corporeal active medical device (e.g. insulin pump, defibrillator, etc.) have been cautioned that some types of electronic equipment might interfere with the operation of such devices. It is recommended that the handpiece and cables be kept at least 6 in (15 cm) away from any pacemaker or pacemaker leads during use. Clinicians should contact the pacemaker manufacturer or the patient's primary physician for detailed information about the pacemaker.

 CAUTION!	• During boil-water advisories, this product should not be operated as an open water system (i.e. connected to a public water system). A Dental Healthcare Professional should disconnect the system from the central water source. When the advisory is cancelled, flush all incoming waterlines from the public water system (e.g. faucets, waterlines and dental equipment) in accordance with the manufacturer's instructions for a minimum of 5 minutes. • Use of the BioSonic US200 Ultrasonic Scaler may irritate demineralized tooth structure, hypersensitive areas, veneers, cast crowns and implants. • Use of the BioSonic US200 Ultrasonic Scaler is only intended to function as intended by the manufacturer. Observe the following to avoid damaging the device or impairing performance: • Do not restrict air flow. • Provide adequate ventilation. • Do not subject unit to shock or impact. • Do not place the unit on or next to a radiator or other heat source. • Do not attempt to sterilize the handpiece cord and core assembly. Only the handpiece sheath may be sterilized. • Do not use petroleum based solvents, iodophors or phenolic based products. (iodophors and phenolics can stain the surface of the unit.) Clean up all liquid spills immediately.
--	---

VII. Order Information

Unit

Catalog No.	Input Voltage/Current	Frequency	Power	Color
60034537	120V/1.5A	50/60 Hz	60W	Dark Gray
60043042	120V/1.5A	50/60 Hz	60W	Purple
60043043	120V/1.5A	50/60 Hz	60W	Cyan
60043044	120V/1.5A	50/60 Hz	60W	Dark Teal
60043045	120V/1.5A	50/60 Hz	60W	Pink

Inserts: Frequency 25 kHz

Catalog No.	Description	Qty
60034383	#10 Universal	1
60034373	#10 SuperSoft	1
60034379	#10 OptiMist	1
60034385	#10 Universal Slim	1
60034377	#10 Slim SuperSoft	1
60034387	#10 Slim OptiMist	1
60034375	#10 Slim SuperSoft (Optimist)	1
60034381	#1000 Universal Triple Bend	1

Inserts: Frequency 30 kHz

Catalog No.	Description	Qty
60034384	#10 Universal	1
60034374	#10 SuperSoft	1
60034380	#10 OptiMist	1
60034386	#10 Universal Slim	1
60034378	#10 Slim SuperSoft	1
60034388	#10 Slim OptiMist	1
60034376	#10 Slim SuperSoft (Optimist)	1
60034382	#1000 Universal Triple Bend	1

Replacement Parts:

60036143	Replacement O-Ring Gaskets	12
60036144	Water Filter Replacement Kit	10
60034538	US200 Replacement Core O-Ring	3
60034539	US200 Handpiece Sheath	1
60034540	US200 External Power Supply	1
60034542	US200 Footswitch	1
60034543	US200 Handpiece Assy	1
60036145	US200 Water Line	1

VIII. Troubleshooting

Generally, check all lines and connections to and from the unit. Check the settings on the unit's knobs. Although service and repair of the BioSonic US200 ultrasonic scaler should be performed by authorized personnel, the following are some basic troubleshooting procedures that may help avoid unnecessary service calls.

A. Unit will not operate:

LED Status	Potential Solution
Off	Ensure the power cord is fully seated in the receptacle on the power adapter.
	Ensure the power adapter cord is fully seated in the DC receptacle on the back of the unit.
	Ensure the power cord is plugged into an active and appropriate AC receptacle.
	Ensure the power knob is in the ON position by checking that the power knob is not rotated fully counterclockwise. Knob should "click" when rotated clockwise, turning power on.
	Ensure the power adapter being used is the one shipped with the unit.
Green	Release the footswitch and reengage. The unit will not activate if the footswitch is depressed when the unit is turned on.
	Ensure the handpiece connector is fully seated in the receptacle on the front of the unit. A "click" will be felt when it is properly inserted.
Blinking Green and Orange	Ensure the footswitch is fully inserted in the receptacle in the back of the unit. A "click" will be felt when it is properly inserted. Inspect the footswitch and cord for any obvious damage.
	Ensure an appropriate 25 kHz or 30 kHz insert is fully seated in the handpiece.
	Ensure adequate water flow to the handpiece.
	Replace the insert. It may be worn, clogged, or damaged.

B. Unit will operate but LED is not solid green or orange:

LED Status	Potential Solution
Blinking Green (in Normal Mode) or Blinking Orange (in Turbo Mode)	Ensure adequate water flow to the handpiece. Ensure water supply to the back of the unit is properly connected and not kinked or pinched. Ensure water filter is not clogged. If clogged, replace with new filter. Adjust water control knob on front of unit until an appropriate mist is observed.
	Ensure an appropriate 25 kHz or 30 kHz insert is fully seated in the handpiece.
	Reduce power setting.
	Replace the insert. It may be worn, clogged, or damaged.

Blinking Green and Orange	Release the footswitch and reengage.
	Ensure the handpiece connector is fully seated in the receptacle on the front of the unit. A "click" will be felt when it is properly inserted.
	Rotate power knob fully counter-clockwise to turn power off, wait 10 seconds, and rotate power knob clockwise to turn power back on.
	Replace the insert. It may be worn, clogged, or damaged.

C. No water at insert:

1. Remove the insert and ensure that the water control is in the fully open position.
2. Check that the main shut-off valve is open.
3. Check the disposable in-line water filter for possible clogging.
4. Check the recommended customer supplied in-line water filter for possible clogging.
5. Make sure that the water lines are not kinked.
6. Make sure that the water knob is turned to the "open" position (clockwise).
7. Check the insert tip for possible clogging

D. Unit sprays water, but the insert does not vibrate:

1. Try a different insert. Sometimes a worn out and/or damaged insert is not detectable to the naked eye.

E. Unit stops operating:

1. If an over-current or over temperature condition is detected, power and waterflow to the handpiece is cut off and the LED will blink green and orange. If this occurs, release and then re-engage the footswitch.

IX. Specifications**A. Unit**

Catalog No.	Input Voltage/Current	Frequency	Output Power	Color
60034537	120V/1.5A	50/60 Hz	60W	Dark Gray
60043042	120V/1.5A	50/60 Hz	60W	Purple
60043043	120V/1.5A	50/60 Hz	60W	Cyan
60043044	120V/1.5A	50/60 Hz	60W	Dark Teal
60043045	120V/1.5A	50/60 Hz	60W	Pink

Operating Frequency: 25/30 kHz

Overall Dimensions: 5 x 7 x 4.5 inches
(12.7 x 17.8 x 11.4 cm)

Weight: 2 lbs (0.9 kg)

B. Inserts **C. Explanation of Symbols**

	Caution/Warning
	Type B Applied Part

	Footswitch
	Water Connection
	Refer to Instruction Manual
	Burn Hazard
	Electric Shock Hazard
	Temperature Limitation
	Humidity Limitation
	Atmospheric Pressure Limitation
	Consult the instructions for use
	Medical Device Indicator: this product meets the definition of a medical device
	Sale Restricted: The purchase of the device is restricted to licensed healthcare practitioners only
	Date of manufacture
	Serial Number
	COLTENE Catalog Number
	The entity that manufactures the medical device
IPX2	Protection Class – IPX2 (Unit) IPX2 Classification of ingress of water
IPX8	Protection Class- IPX8 (footswitch only) IPX8 Classification of ingress of water

X. Classifications

Protection against electrical shock:

Base Unit - Class II

Inserts - Type B

Protection against ingress of water: IPX8 (footswitch) / IPX2 (unit)

Mode of operation: Continuous

Equipment not suitable for use in the presence of flammable mixtures.

XI. Disposal of unit/waste products

In accordance with local and state laws.

XII. Limited Warranty

Dent4You AG products are carefully manufactured to meet stringent quality assurance requirements. Dent4You AG products are manufactured from new parts or new and serviceable used parts. Regardless, Dent4You AG's warranty terms apply. This product has been developed specifically for use in dentistry and is intended to be operated only by qualified dental professionals in accordance with the instructions contained in this guide. However, notwithstanding anything contained herein to the contrary, the user shall at all times be solely responsible for determining the suitability of the product for the intended purpose and the method of its use. Any guidance on application technology offered by or on behalf of the manufacturer, whether written, verbal or by demonstration, shall not relieve the dental professional of his/ her obligation to control the product and to make all professional judgments regarding its use.

Dent4You AG warrants this product to be free from defects in workmanship and materials for a period of twenty-four (24) months from the date of original purchase. During such warranty period, Dent4You AG undertakes to repair or replace this product, at its option, if such product does not function properly under normal use and service, which malfunction is due solely to a defect in workmanship or material.

Any claim for damage or breakage to the product in transit should be made to the carrier promptly upon discovery. Dent4You AG does not warrant the product against shipping damage.

This warranty does not cover any of the following:

- Damage caused by accident, abuse, misuse, neglect, alteration, transit, failure to follow manufacturer's instructions, lack of reasonable and proper maintenance and care, or Acts of God.
- Normal wear and tear (i.e.: bulbs, fuses, etc.).
- Minor imperfections that do not materially affect the product's performance.
- Use of the product by anyone other than a licensed dental practitioner or qualified dental hygienist

This warranty shall become null and void if any repair or servicing of the product is made or attempted by any person not authorized by Dent4You AG, or if any part not authorized by Dent4You AG is used in any repair or servicing.

This warranty is in lieu of all other warranties relating to the product, expressed or implied, including, without limitation implied warranties of merchantability and fitness for a particular purpose. In no event shall Dent4You AG be liable for direct, indirect tort, exemplary, punitive, special, incidental, or consequential damages arising from the possession or use of the product, including but not limited to loss of use of the product, lost revenues or profits, or the cost of substitute equipment. The limitations provided for herein shall inure to the benefits of Dent4You AG's suppliers and subcontractors.

This Certificate of Limited Warranty is valid only if it has been completed with date of purchase and dealer's name.

If this product becomes defective during the warranty period, it must be delivered together with this Certificate of Limited Warranty and a written warranty claim to the Dent4You AG dealer from which it was originally purchased. All warranty claims must include a detailed written description of the claimed defect. The purchaser will be responsible for all costs and risks of transit, and we therefore suggest that the product be carefully packed and that it be shipped to the dealer insured for its full value.

Register the warranty at warranty.coltene.com.

XIII. EMC Information, Warnings and Considerations:

1. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instruction manual, it may cause interference to radio communications. Operation of this equipment in a manner inconsistent with this manual is likely to cause interference in which case the user, at their own expense, will need to take whatever measures may be required to correct the interference.
2. Persons fitted with cardiac pacemakers have been cautioned that some types of electronic equipment might interfere with the operation of a pacemaker. Although no instance of pacemaker interference relating to a BioSonic Ultrasonic Scaler has ever been reported to Dent4You AG Inc., we recommend that the handpiece and cables be kept at least 6 inches (15 cm) away from any pacemaker and pacemaker leads during use. If use of this product is in question, consult patient's cardiologist.
3. Use of BioSonic US200 cables and accessories, other than those provided by Dent4You AG, may result in increased emissions or decreased immunity.
4. Listed cables:
 - Handle & Cord Assembly – Part # 60034543
 - Footswitch & Cord Assembly – Part # 60034542

5. Electromagnetic Compatibility:

The following are guidance and manufacturer's declarations regarding electromagnetic compatibility for the BioSonic US200.

Conforms to UL STDs:	Conforms to CSA STDs:
60601-1	C22.2 No. 60601-1
60601-1-6	C22.2 No. 60601-1-6
60601-2-60	C22.2 No. 60601-2-60

NOTE The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment

5.1 EN/IEC 60601 - 1 -2 Table 1

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The BioSonic US200 is intended for use in the electromagnetic environment specified below. The customer or the user of the BioSonic US200 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - Guidance
RF emissions CISPR 11	Group 1	The BioSonic US200 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The BioSonic US200 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

5.2 EN/IEC 60601 - 1 -2 Table 2

Guidance and Manufacturer's Declaration – Immunity			
The BioSonic US200 is intended for use in the electromagnetic environment specified below. The customer or user of the BioSonic US200 should ensure that it is used in such an environment.			
Immunity Test	IEC60601 test level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	+2kV for power supply lines +1kV for input/output lines 100 kHz repetition frequency	+2kV for power supply lines +1kV for input/output lines 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV Line-to-line ±0.5 kV, ±1 kV, ±2 kV Line-to-Ground	±0.5 kV, ±1 kV Line-to-line ±0.5 kV, ±1 kV, ±2 kV Line-to-Ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BioSonic US200 requires continued operation during power mains interruptions, it is recommended that the BioSonic US200 be powered by an uninterruptible power supply or battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30A/m 50/60 Hz	30A/m 50/60 Hz	Power frequency magnetic fields should be that of a typical commercial or hospital environment

5.3 EN/IEC 60601 - 1 -2 Table 2




Guidance and Manufacturer's Declaration – Immunity			
The BioSonic US200 is intended for use in the electromagnetic environment specified below. The customer or user of the BioSonic US200 should ensure that it is used in such an environment.			
Immunity Test	IEC60601 test level	Compliance Level	Intended Electromagnetic Environment
Conducted RF IEC 61000-4-6	3 Vrms 6Vrms (In ISM Bands) 150 kHz to 80 MHz 3V/m	3Vrms 6Vrms (In ISM Bands) 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the BioSonic US200 including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ $d = 2,3\sqrt{P}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitter as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	80MHz to 2.7GHz	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	

Proximity fields from RF wireless communications equipment

Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Moducation ^{b)}	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse moducation ^{b)} 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1,720	1,700 - 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 - 2,570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5,240	5,100 - 5,800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
5,500						
5,785						

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©2024
Coltene/Whaledent, Inc.