

NEW COLORS

BioSonic® US200

Bold colors. Bright smiles.







Experience the ultimate in ultrasonic scaling

COLTENE's BioSonic US200 ultrasonic scaler delivers the power and reliability today's hygienists need to clean deeply, maximize chairside efficiency, and provide a comfortable treatment experience for patients. And now it comes in four head-turning new colors.

AVAILABLE COLORS








KEY FEATURES

-  Increased power range with smooth transitions
-  Easy-to-see LED indicator light
-  Precise water flow control
-  Lightweight, ergonomic design
-  Automatic insert frequency detection
-  Two-stage foot pedal for easy power adjustments



BIOSONIC SCALER INSERTS

	UNIVERSAL #10			TRIPLE BEND #1000	
	25 kHz	60034383		25 kHz	60034381
30 kHz	60034384		UNIVERSAL #10 SUPERSOFT		
25 kHz	60034379		25 kHz	60034373	
30 kHz	60034380		SLIM #10 SUPERSOFT		
25 kHz	60034385		25 kHz	60034377	
30 kHz	60034386		SLIM #10 OPTIMIST SUPERSOFT		
25 kHz	60034387		25 kHz	60034375	
30 kHz	60034388		30 kHz	60034376	

Colors that pop. Performance that delivers.



ORDERING INFO

WHAT'S INCLUDED	TECHNICAL INFO	COLOR	PART NO.
BioSonic US200 Ultrasonic Scaler Unit Footswitch Scaler Handpiece Assembly Universal Power Supply Water Supply Tubing Ultrasonic Scaler Insert Sterilizable Handpiece Sheaths, 3 pcs Handpiece Sheath O-Rings, 3 pcs Water Line Filters, 3 pcs	Overall Dimensions: 5 x 7 x 4.5 in (12.7 x 17.8 x 11.4 cm) Weight: 2 lbs (0.9 kg) DC Power Supply: 24 VDC - 60W min	BLACK	60034537
		PURPLE	60043042
		CYAN	60043043
		DARK TEAL	60043044
		PINK	60043045

**SCAN TO LEARN MORE
& SEE SPECIAL OFFERS**

COLTENEGROUP Questions? Get in touch with your local COLTENE Territory Manager today.

COLTENE
coltene.com

