

BioSonic® Solution Concentrates

The solution for your cleaning needs

For Use With All Ultrasonic Cleaning Units

- ✓ Super concentrated
- ✓ Super effective
- ✓ Super convenient



UC30 Super Concentrate General Purpose Ultrasonic Cleaning Solution

- Dissolves and removes blood, saliva and microscopic particles
- Perfect for cleaning all dental instruments – even plastic and rubber
- 16 oz. bottle yields 8 gallons of solution



UC32 Enzymatic Ultrasonic Cleaning Concentrate

- Minimizes hand scrubbing by removing protein soils (dried blood, tissue, saliva)
- Effectively breaks down and loosens hardened debris
- 8 oz. bottle yields 8 gallons of solution



UC34 Plaster and Stone Remover Powder Concentrate

- Concentrated powder for easy cleaning of attachments, articulators, instruments, castings and items that cannot be sandblasted
- 1.85 pound jar yields 3 gallons of solution



UC38/42 Germicidal Ultrasonic Cleaning Concentrate

- Kills HIV virus and tuberculosis in the bath while it cleans instruments
- EPA registered
- Excellent as a pre-soak
- UC38 = 8 oz. bottle yields 8 gallons of solution
- UC42 = 16 oz. bottle yields 16 gallons of solution



UC39 Ultrasonic Cement Remover Concentrate

- Removes virtually all types of temporary and permanent cements from instruments and prostheses
- 16 oz. bottle yields 32 applications



Exclusive MeterDose® Bottle

- Dispenses exactly the amount needed; No waste, no spills, no mess
- Simply squeeze, tilt and pour
- No bulky gallon jugs to store; save up to 85% in storage space
- Fast acting; no need to wait for tablets or powder to dissolve

Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, OH 44223 / USA
Tel. USA & Canada 1 800 221 3046
+1 330 916 8800
Fax +1 330 645 8704
info.us@coltene.com

Call Your Local Dealer Today
& Ask For BioSonic Solutions

www.coltene.com
COLTENE

P/N 30002968C